

Vaccine Temperature Log (Fridge and Freezer)

Clinic

Fridge (make / ID)

Freezer (make / ID)

Month and year

CDC: fridge 2 to 8 °C (36 to 46 °F); freezer -50 to -15 °C (-58 to +5 °F). Check and record min/max at the start of each workday.

Day	Time	Fridge now	Fridge min	Fridge max	Freezer now	Freezer min	Freezer max	Initials	Excursion? Action taken
1									
2									
3									
4									
5									
6									
7									
8									
9									
10									
11									
12									
13									
14									
15									
16									
17									
18									
19									
20									
21									
22									
23									
24									
25									
26									
27									
28									
29									
30									
31									

Reviewed by: _____ Date: _____ Unit (circle): °F / °C